

Parts Catalogue

Therapy Products




Parts

If you need further information regarding the parts for your Laerdal products, please contact us. Trained personnel are available to answer your questions and provide additional technical data.

All parts listed for a specific product might not be available in all markets

Product specifications are subject to change without notice. Actual colours on the parts might differ from the colours in this catalogue.

For online purchase please visit
www.laerdal.com



Contents

HeartStart MRx 4-7
HeartStart FR – FR2+ 8
HeartStart HS1-FRx 9
HeartStart FR3 10
Laerdal Silicone Resuscitators 12-13
Pocket Mask 14
Paediatric Pocket Masks 14
Stifneck Select 14
SpeedBlocks Head Immobilizer 14
Baxstrap Spineboard 15
V-VAC 15
Laerdal Suction Unit 16-17
Laerdal Compact Suction Unit LCSU4 18

Contents



M3541A
Carry case for HeartStart MRx, incl. 3 accessory pouches/shoulder strap



989806100851
MRx Ambulance Wall Mount



M3538A
MRx battery



- 40457C** Recording Paper Rolls For Cm100 (Pack Of 10 Rolls)
- 40457D** 50mm printing paper (HS 3000QR, HS 4000, MRx) (80 rolls)
- 989803138171** MRx 75mm printer paper 10 rolls
- 989803138181** MRx 75mm printer paper 80 rolls



948003 Power cable



948002 DC Power Supply for MRx






- M3725A** Test load for use with M3508A cable (Plug connector)
- M1781A** Test load for use with M3507A cable (Barrel connector)
- 94 50 02** HeartStart FR to Symbio Tester adapter cable (not for patient use)


Temperature probes for HeartStart MRx (M3536A/M3535A):

- 21078A** Reusable Skin Surface Temperature Probe
- 21082A** Reusable connector cable for disposable temperature probes.
- 21090A** Disposable Oesophageal/Rectal Temperature Probe
1 box = 20 probes. (Requires 21082A)
- 21091A** Disposable Skin Surface Temperature Probe
1 box = 20 probes. (Requires 21082A)





ECG Cables

ECG CABLES	
	<p>989803147691 Trunk Cable (1,3 meters) (IEC) (standard from Q3/06)</p> <p><i>The old long trunk cable (2,7m) can be ordered as sparepart (M1949A)</i></p>
	<p>M3526A 3-lead ECG patient cable (snap) (1,6 meter) (IEC)</p>
	<p>M3529A 7-lead ECG patient cable (snap) (IEC) (To be used together with M3526A)</p>

ECG 5 LEAD ECG	
	<p>M1971A 5-lead ECG patient cable (extremities. + V1)</p> <p>M1978A 5-lead ECG patient cable – chest (V2 to V6) (To be used together with M1971A)</p>

SpO₂

<p>M1191A SpO₂ Adult (>50 kg)</p> 	<p>M1194A SpO₂ Ear clip (>50 kg)</p> 
<p>M1192A SpO₂ Child (15-50 kg)</p> 	<p>M1195A SpO₂ Infant (4-15 kg)</p> 
<p>M1941A SpO₂ Adapter extension cable (2 meter)m</p> 	

Nelcor Disposable SpO ₂ Sensor	
	<p>M1943A SpO₂ adapter to Nelcor single use probe (M1902/4B) (1,1 meter)</p>
	<p>M1904B Adult (>50kg) M1903B Pediatric (15-50kg) M1902B infant (4-15kg)</p>

NBP Antimicrobial Reusable Cuffs



- M4552B Easy Care Cuff, Infant (10 - 15 cm Circ.)
- M4553B Easy Care Cuff, Pediatric (14 – 21.5 cm Circ.)
- M4554B Easy Care Cuff, Sm Adult (20.5 - 28.5 cm Circ.)
- M4555B Easy Care Cuff, Adult (27.5 – 36 cm Circ.)
- M4556B Easy Care Cuff, Adult XL (27.5 – 36 cm Circ.)
- M4557B Easy Care Cuff, Large Adult (35 – 45 cm Circ.)
Cuff cover length 66.5 cm.
- M4558B Easy Care Cuff, Large Adult XL (35 – 45 cm Circ.)
Cuff cover length 83 cm



- M1598B
1,5 meter NIBP tubing
- M1599B
3 meter NIBP tubing

Reusable Comfort Cuffs



- M1572A Pediatric size cuff
- M1573A Adult small size cuff
- M1574A Adult size cuff
- M1575A Large Adult size cuff
- M1576A Thigh size cuff
- M1578A Multi-Patient Reusable Cuff Kit
(M1573A, M1574A, M1575A, M1576A)


HeartStart MRx Cat.No. M3536A

External Defibrillation Paddles



M3543A
Water resistant external anterior/ anterior paddle set with patient contact indicator (PCI)

Pads Cable



M3508A
Pads Cable – (Plug connector)



M3507A
Pads Cable – (Barrel connector)



M4763A
Therapy/CPR Cable

Adapters

05-10200
Pads Adapter
(Barrel connector to plug style)

05-10100
Zoll™ Pads Adapter
Can be used with Zoll 1200, 1400, 1600 and M-Series



05-10000
QUIK-COMBO Pads Adapter.
Can be used with Medtronic Physio-Control LifePak 9, 10C, 11, 12.



989803146981
256MB MRx Data
Card and Tray



M3854A
Datacard and tray



Kat. No	Description	CO ₂ Only	O ₂ /CO ₂	Intubated/ LMA or Combitube	Nasal Only	Oral/Nasal	Duration for use	Primary Applications
M2520A	Pediatric Smart CapnoLine™ O ₂		.			.	8H	10-20 kg
M2522A	Adult Smart CapnoLine™ Plus O ₂		.			.	8H	>55 kg
M2524A	ETCo ₂ Pediatric Smart Capnoline Single purpose	.				.	8H	10-20 kg
M2526A	ETCo ₂ Adult Smart Capnoline Single purpose	.				.	8H	>55 kg
M1920A	Adult/Pediatric FilterLine® Set	.		.			6H	>2 kg
M1923A	Infant/Neonatal	.		.			72H	>2 kg

Dual purpose = a probe that enable monitoring of CO₂, and enable supply of O₂.
Humidified = can be used up to 24 hours. All others can be used to 8 hours.



M2520A
Peadiatric Smart Capnoline O₂



M2522A
Adult Smart Capnoline



M1920A
Adult/Pediatric MicroStream
Filterline


M2522A
Adult Smat Capnoline O₂

Therapy (Shock & Pacemaker) electrodes

Electrodes for FR2, MRx, HS 4000



989803158211 Adlt Defibrillation Pads (1 Pair)
 989803158221 Adult Defibrillation Pads (5 Pairs)
 M3713A
 Adult/Child Plus with plug style connector



M3717A
 Infant Plus with plug style connector

Electrodes for MRx, HS 4000
 (with barrel style connector)



M3502A Adult Multi Elect For Use With
 Cm100/Heartstart4000 (Pack Of 10)



9851-008
 Symbio CS1201 (12-lead rhythm simulator tester)

9801-008
 Symbio CS301 (3-lead rhythm simulator tester)



989803162401
 MRx Replace CPR Meter



801-10850 Patient
 Adhesives (10x3)

NOT ILLUSTRATED

989806100852 Lower Mount "A" Bracket	989803135341 4 Bay Battery support system for 2 sealed lead acid and 2 lithium ion batteries
989806100853 Lower Mount "B" Bracket	861440 HeartStart Telemedcare System Classic Care Edition (English Only)
989806100854 Upper Mount, Top Bracket	989803146981 256 MB MRx Data card and Tray
989803163291 Patient Adhesive (10pc.) DP2 Defibrillation electrodes (2 pair)	989803158081 5-lead set, limb, IEC
M1191B Reusable Adult SpO2 finger sensor	989803158091 5-lead set, chest, IEC
M3716A Adult/Child Radiolucent with plug style connector	989803158661 Repl. Pads/CPR Cable
M3718A Adult/Child Radiotransparent (Reduced Skin Irritation) with plug style connector	989803172501 MRx Black Soft Carry Bag
M3719A Infant Radiotransparent (Reduced Skin Irritation) with plug style connector	989803176161 5 lead set, snap, shield AAMI, Chest
21075A Esophageal/rectal temp probe, Adult Reusable	989803176171 5 lead set, snap, shield AAMI, Limb
M3527A 7-lead chest lead set, snap, AAMI 321040	989803176181 5 lead limb lead wires
989803135301 2-BAY ANALYZER FOR HS Li-ION BATTERIES	989803176191 5 lead chest lead wires
989803135331 Battery charger: 4 Bay Battery support system for 4 lithium ion batteries	M1193A Neonatal Hand/Foot SPO2 Sensor
	M1530A 5 Lead ECG Cable, IEC
	M1663A 10-ELECTRODE CABLE SETS - 5 WIRE
	M2202A Adult radiotranslucent foam solid gel electrode Size:2"(50mm) diameter
	M3501A Adult Pads with barrel connector
	M3504A Infant Pads with barrel connector
	M3539A MRx Power Module
	M3542A External Paddles
	M3545A MRx Data Card and Tray
	M5521A MRx handle (Green)
	M5522A MRx handle (Blue)
	M5525A MRx Handle (grey)
	M5529A Heartstart MRx DC Power Module
	M1510A 3 lead ECG Patient Trunk Cable, IEC
	M1615A Shielded 3-Lead Set, Snaps, IEC
	M1633A Shielded 5-Lead Set, Grabbers, IEC
	M1635A Shielded 5-Lead Set WITH Safety Conn.
	M1669A 3 LEAD TRUNK CABLE AAMI/IEC 2.7m
	M3528A 3-electrode limb lead set, snap, IEC
	M3549A Wide Bed Rail Hook (Metal)
	M4737A Keyboard Overlay MRx
	M5530A Combiner plug for 3-wire lead set
	989803166021 Adult / Child Pre-connect
	989803158061 5-lead set, limb, AAMI
	989803158071 5-lead set, chest, AAMI
	M4612A Adult solid get ECG electrode Shape: teardrop
	M4613A Adult solid get ECG electrode Shape: teardrop
	M4761A Q-CPR Compression Sensor
	M4762A CPR Assistent Replacement Pads (10)
	989803185181 Black Soft Carry Case
	989803176541 HeartStart MRx Quick Disconnect DC Power Cable

HeartStart FR - FR2+ Cat. Nos. M3841A, M3840A

Batteries



M3863A
Long-Life LiMnO2
Battery. Disposable
and recyclable. For
FR2 and FR2+ only.



M3854A
Datacard and tray



M3870A Infant/Child
Reduced-Energy
Defibrillator Pads



07-10900
Training pads



40457C
Recording Paper Rolls
(Pack Of 10 Rolls)



989803158211 Adult Pads (1 set)

989803158221 Adult Pads
(1 case = 5 sets)

M3713A
Multifunction Electrode Pads
Adult Plus (1 case = 10 sets)



M3864A
Training & Administration Pack
for FR2 series Defibrillator



05-10100
Zoll™ Pads Adapter
Can be used with Zoll 1200, 1400,
1600 and M-Series



05-10000
QUIK-COMBO Pads Adapter.
Can be used with Medtronic Physio-
Control LifePak 9, 10C, 11, 12.

NOT ILLUSTRATED

945021	Label status AED Trainer 2 (pkg. 10)
945025	Board AED Trainer 2
945026	Speaker with cable
945027	Cable AED Trainer 2
945055	Programming kit for AED Trainer 2 (optional software and cable)
945060	Technical manual
DP2	Defibrillation electrodes (2 pair)
M3871A	FR2 Paediatric Training Pads
68-PCHAT	Fast Response Kit
941308	Softpack, FR2 Army Green
941320	FR2 Rescue Rucksack
941404	Compact Flash Card Reader USB
M3853A	Data Card Tray only
M3869A	Vinyl Carrying Case for FR2 AED



941311
Carry Case (Hard Pelican)

Adapters

05-10200
Pads Adapter
(Barrel connector to plug style)

Discontinued model

HeartStart HS1 Cat. Nos. M5066A, M5067A


Shock Electrodes



M5071A
Adult SMART Pads
Cartridge (1 Set)



M5072A
Infant/Child SMART
PADS Cartridge (1 set)



M5073A
Adult Training Pads
with Replaceable Seal

Battery

M5070A
(For HS1 and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



NOT ILLUSTRATED

- M5075A Infant/child training pads cartridge
- M5093A Replacement adult training pads
- M5094A Replacement infant/child training pads
- M5089A External manikin adapters (10 on a roll)
- 941320 Rescue Rucksack
- 941343 HeartStart Wall Box with Alarm
- 941345 Heartstart kit Barrier, wipe & gloves
- 989803139531 Carry case
- 989803143051 Case Capture palm SW
- 989803170911 AED Green Awareness Placard
- 989803170931 AED Wall Sign, Green
- M3859A Secure Pull Seal Wall Mount
- M5075A HeartStart HS1 Carry Case
- M5076A HeartStart HS1 Slim Carry Case
- M5090A Adult Pad Placement Guide

HeartStart FRx Cat. Nos. M861304, 861305

Carry Case

941360
FRx carrying case made of wipeable urethane.



Battery

M5070A
(For HS1 and FRx) This battery is a disposable, lithium manganese dioxide, long-life battery that will typically last for four years with the defibrillator in standby mode.



NOT ILLUSTRATED

- 940314 Bracket without Lock
- 941345 Heartstart Kit Barrier, wipe & gloves
- 989803139251 HeartStart FRx Carry Case
- 989803139281 HeartStart Pad Placement Guide, Infant / Child
- 989803170891 Wall Mount Bracket
- 989803170911 AED Awareness Placard(Green)
- 989803170931 AED Wall Sign, Green
- M3858A AED Wall sign
- M3859A Secure Pull Seal Wall Mount 10 pack
- M5090A Adult Pad Placement Guide

Shock Electrodes



989803139261
SMART Pads II




989803139311
Infant/child key used in infant/child defibrillation with the HeartStart FRx

989803139271
Training Pads II (1 pair)

989803139291
Replacement Training Pads TL (1 pair)





989803150191 FR3 Training Pack



945006 Cable AED to Manikin FR/FR2/FR3



989803149981 SMART Pads III, 1 set
989803149991 SMART Pads III, 5 set



801-10850 Patient Adhesives (10 x 3)



989803150111 Fast Response Kit



989803149941 Q-CPR Upgrade

NOT ILLUSTRATED

- 989803150181 Replacement Pads HS Training Pads III
- 989803150201 Interconnect Cable HS Training Pads III
- 989803149951 Cable link,FR3 to CPR meter
- 989803150041 3-Lead ECG cable,FR3 AAMI
- 989803150051 3-Lead ECG Cable,FR3 IEC
- 989803150241 FR3 Rechargeable Battery For Clinical Use
- 989803157531 Laerdal System Case, Rigid
- 989803162231 CPR meter cradle
- 989803179171 Laerdal System Case Soft, FR3 (W/O Auto-on)
- 989803179191 Laerdal small Soft Case FR3 (W/O Auto-on)
- 989803184831 Laerdal small Soft Case FR3 (W/O Auto-on)
- 989803150241 FR3 Rechargeable Battery For Clinical Use

Laerdal Silicone Resuscitators Cat. Nos. 870050/51/52/53/55,

560200 Non-Rebreathing Valve Iv



540105
Disk membranes
(pkg. 10)

540103
Lip valve



870150
Adult Bag 1600MI

851250 Patient valve w/press. reg.



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252
Pressure Relief Valve 35Cm H2O



860150
Paediatric Ventilation bag,
500 ml

851350 Patient valve w/press. reg. and lock clip



540105
Disk membranes
(pkg. 10)

540103
Lip valve

851252
Pressure Relief Valve 35Cm H2O

851103
Lock clips
(pkg. 10)



850150
Preterm Bag 240MI



530400 Airways, set of 4

TUBE AND CONNECTORS



871000 Silicone extension
tube (28 cm)



OD15 mm
OD6 mm
ID15 mm
OD22 mm
850900
Manometer
connector



OD30 mm
850500
Expiration
diverter



875400
Intake/Reservoir
Valve



551901 O₂ reservoir,
0.6 litre



531901 O₂ reservoir,
2.6 litres

NOT ILLUSTRATED	
54010733	High Efficiency Filter for LSR/TheBag(IE)(2.5pk)
531906	O2 Reservoir 2.6Ltr (Pack Of 50)
551906	O2 Reservoir 0.6Ltr (Pack Of 50)
85015005	Preterm ventilation bag, 240ml (J)
86015005	Paediatric ventilation bag, 500ml (J)
87015005	Adult ventilation bag, 1600ml (J)
87540000	Intake Reservoir Valve (Nordic)
87540033	Intake Reservoir Valve (IE)
87540040	Intake Reservoir Valve (EUR)



860300 Display Case Paediatric
870600 Display Case Complete Ad.Mod.
850700 Display Case Preterm



875400 Intake/Reservoir Valve



510103
Cap for LSR Intake
Valve (pkg. 3)



510404
Intake
membranes
(pkg. 10)



871950 Flap valves
(pkg. 2)



531907
Intake Valve Adapter
(23 mm OD)



870120
Hanging loop



511700
Wall bracket



871300
Attachment ring for
standard head strap



870400 Head strap
w/attachm. ring

Pocket Mask Cat. No. Paediatric Pocket Mask Cat. No.

NOT ILLUSTRATED

820410	One Way Valve Single
820411	Valve Pocket Mask 100/pkg.
820700	Update kit Pocket Mask
830110	Head Strap P/Mask (Pack Of 10)
820610	Pocket Mask Mask Filter Single
820611	Pocket Mask Mask Filter 10 Pack

Stifneck Select Collars Adult & Cat. No. Stifneck Extrication Collars Cat. No.

NOT ILLUSTRATED

980700	Stifneck Carrying Bag
--------	-----------------------

SpeedBlocks Head Immobilizer Cat. No.



982100 PadPack
Alignment Pads (Pack 10)

NOT ILLUSTRATED

983060	SpeedBlocks Universal Base (Qty.1)
983092	SpeedBlocks Block Set (Qty.1)
983096	SpeedBlocks Strap and Pad Replacement Set

BaxStrap Spineboard Cat. No.



982100 PadPack
Alignment Pads (Pack 10)



982512 Premium Strap
w/sewn loop ends



982513
Premium Strap w/
speed clip ends



982517 Spider Straps
for BaxStrap

NOT ILLUSTRATED

- | | |
|-----------|---|
| 701-00040 | Laerdal Pedi-Pad
Backboard Pad (Int) 6pk |
| 980800 | Fullcover Carry Bag
Baxstrap Spineboard |

V-VAC Cat. No.



NOT ILLUSTRATED

- | | |
|--------|---------------------------------|
| 886010 | 300ml Disposable |
| 985001 | V-Vac Replace Cartridge |
| 985002 | Adaptor Tips
(Pack Of 4) |
| 985003 | Male Connectors
(Pack Of 10) |
| 985004 | Catheters (Pack Of 4) |
| 985101 | V-Vac Training Cartridge |
| 985200 | V-Vac Carry Bag
(Yellow) |
| 985030 | V-VAC Handle
Spare Part |

Laerdal Suction Unit Cat. Nos. 780000, 780020, 780030

Reusable Canister System Cat. No. 780000

650113
Suction catheter
adapter
(pkg. 10)



770410
Suction tube
150 cm
w/o tip



781200
Aerosol Filter for LSU
Reusable Canister

793500
Water
canister



780435
Holder for
water canister

784000 Reusable canister


781006
Angled
connectors
(pkg. 10)



781002
Float Ball
(pkg. 10)



784007
Vacuum
seal
(pkg. 10)



784008
Gaskets
(pkg.10)



Serres Canister System Cat. No. 780030



57151 Serres Suction Bag 1000ml
57151-2 Serres Suction Bag (Pk Of 36)



57300
Serres Canister 1000. ml



781206
Serres Vacuum Connector



5833181 Serres Suction Tubing

Serres Canister System Pre 2014



781203 Serres High Efficiency
Filtration Kit LSU



781204
Serres Vacuum Tubing



780450
Serres Canister Holder

BEMIS canister system (US Only) Cat. No. 780020



883000
Disposable collection
canister and suction tubing
(12 pk)



780421
Vacuum connector
tubing f/Bemis



780412 Disposable Tubing



78041333 High efficiency
filter kit Bemis

Parts for all units Cat. Nos. 780000, 780010, 780020



780433 Strap for Tubing



780434 Hook for Tubing Strap



780436 Fasten bracket left/right



780800
LSU NiMH Battery

NOT ILLUSTRATED

- 780437 LSU Fastening Bracket
- 780100 Battery Cable
- 780120 Power Cable
- 780130 Carrying handle
- 780140 Main Rotary Switch
- 780150 Gasket
- 780160xx Cabinet cpl.
- 780170 Vacuum Connector
- 780180 Foam Absorbers LSU
- 780190 Battery Lid
- 780240 Power Board
- 780250 MMI Board
- 780270 Battery Contact w/spring
- 780500 Motor & Pump Assembly
- 780510 Pump wear and tear parts
- 780520 Piston-arm exchange tool
- 780530 Umbrella valves
- 780600 LSU mmHg Display
- 780700 LSU KPa display
- 780110 Technical manual


- Reusable Canister System Pre 20xx**
- 781007 Vacuum
 - 780431 Vacuum canister tubing/
filter tubing 30 cm LSU
 - 781004 Aerosol filter (12 pk)
 - 781005 Aerosol filter (50 pk)

*Parts for authorized service personnel only
For other service parts see service manual*

CANISTER HOLDER



780432 Release arm,
canister holder



780420 Bemis Canister Holder



780451
Serres Canister Holder



780430
LSU Reusable Canister Holder

Serres Conversion Kit



780931xx Serres Conversion Kit

Reusable Conversion Kit



780900xx LSU conversion kit
reusable canister

Serres Alteration Kit



780930xx Serres Alteration Kit

780210
AC power cord US



780220
AC power cord EU

780230
AC power cord UK

780200
DC power cord



782600 Wall Bracket
w/DC-Power Cord

782610 Wall Bracket
w/AC Power Cord US

782620 Wall Bracket
w/AC-Power Cord EU

782630 Wall Bracket
w/AC-Power Cord UK

782640 Wall Bracket
w/o Power Cord



791900
Socket adapter



Laerdal Compact Suction Unit LCSU4



886104
800 ml Disposable Canister
w/o tubing (Qty. 6)



886102
800 ml Disposable Canister
w/tubing



886100
300 ml Disposable Canister
w/tubing

NOT ILLUSTRATED

881151	LCSU4, Service Replacement Unit
881152	LCSU4, Service Replacement Unit



886121 Top Cover, w/ Label for control panel

886221 Top Cover, RTCA, w/ Label for control panel.



886105
Patient Tube 1.8 m (6')
Disposable



886106
Vacuum Tube (for 800 ml)



886111
External Battery Charger



886113
Battery, 12V DC NiMH,
Rechargeable



886122 Bottom Cover



884500
12V DC power cord



886128 AC Power-cord UK



886127 AC Power-cord EU



886126 AC Power-cord US



886123
Battery Cover (door)

Compatible Items

Items that can be used with LCSU 2, 3 and 4

Cat. No	Part Name	2nd Part Name	Can be used with:			Version	
			LCSU 4	LCSU 3	LCSU 2	300ml	800ml
886100	300 ml Disp. Canister	w/tubing (Qty.1)	•	•		•	
886102	800 ml Disp. Canister	w/tubing (Qty.1)	•	•	•		•
886104	800 ml Disp. Canister	w/o tubing (Qty.6)	•	•	•		•
886105	Patient Tube 1.8 m (6')	Disposable (Qty.1)	•	•	•	•	•
886106*	Vacuum Tube	(Qty.1)	•	•	•		•
886108	Wire Stand (for 800 ml)		•				•
886111	AC/DC Adapter Charger	w/o Power-cord	•	•	•	•	•
886112	External Battery Charger	w/o Power-cord	•	•		•	•
886113	Battery, 12V DC NiMH	Rechargeable	•	•		•	•
886114	DC Power-cord		•	•	•	•	•
884500	12V DC Power Cord	(1 each)	•	•	•	•	•
881151	LCSU 4, Replacement unit		•			•	•
881152	LCSU 4, RTCA Repl. unit		•			•	•
886121	Top Cover, w/Label for	control panel	•			•	•
886122	Bottom Cover		•			•	•
886123	Battery Cover (door)		•			•	•
886126	AC Power-cord US		•	•	•	•	•
886127	AC Power-cord EU		•	•	•	•	•
886128	AC Power-cord UK		•	•	•	•	•

* Items are not compatible with 1200 ml canister for LCSU 2.



LCSU 4



LCSU 3



LCSU 2



Laerdal

helping save lives